

EEEEEEEEE XX XX AAAAAAA MM MM PPPPPP
EEEEEEEEE XX XX AAAAAAA MM MM PPPPPP
EEEEEEEEE XX XX AAAAAAA MM MM PPPPPP
EE XX XX AA AA MMMM MMMM PP PP LL
EE XX XX AA AA MMMM MMMM PP PP LL
EE XX XX AA AA MMMM MMMM PP PP LL
EE XX XX AA AA MM MM MM PP PP LL
EE XX XX AA AA MM MM MM PP PP LL
EE XX XX AA AA MM MM MM PP PP LL
EE XXXX AA AA MM MM PPPPPP LL
EE XX AA AA MM MM PPPPPP LL
EE XX AA AA MM MM PPPPPP LL
EE XX XX AAAAAAAA MM MM PP LL
EE XX XX AAAAAAAA MM MM PP LL
EE XX XX AAAAAAAA MM MM PP LL
EE XX XX AA AA MM MM PP LL
EE XX XX AA AA MM MM PP LL
EE XX XX AA AA MM MM PP LLLL
EE XX XX AA AA MM MM PP LLLL
EE XX XX AA AA MM MM PP LLLL

FILEID**LABIOCIN

J 9

LL AAAAAAA BBBBBBBB IIIIIII 000000 CCCCCCCC IIIII NN NN
LL AAAAAAA BBBBBBBB IIIIIII 000000 CCCCCCCC IIIII NN NN
LL AA AA BB BB IIII 00 00 CC III NN NN
LL AA AA BB BB IIII 00 00 CC III NN NN
LL AA AA BB BB IIII 00 00 CC III NNNN NN
LL AA AA BB BB IIII 00 00 CC III NNNN NN
LL AA AA BBBBBBBB IIII 00 00 CC III NN NN
LL AA AA BBBBBBBB IIII 00 00 CC III NN NN
LL AAAAAAAAAA BB BB IIII 00 00 CC III NN NNNN
LL AAAAAAAAAA BB BB IIII 00 00 CC III NN NNNN
LL AA AA BB BB IIII 00 00 CC III NN NN
LL AA AA BB BB IIII 00 00 CC III NN NN
LLLLLLLLLL AA AA BBBBBBBB IIIII 000000 CCCCCCCC IIIII NN NN
LLLLLLLLLL AA AA BBBBBBBB IIIII 000000 CCCCCCCC IIIII NN NN

: KW MIST = 1
.TITLE LABIOCIN - LABIO Connect-to-Interrupt Module
.IDENT 'V04-000'

* COPYRIGHT (c) 1978, 1980, 1982, 1984 BY
* DIGITAL EQUIPMENT CORPORATION, MAYNARD, MASSACHUSETTS.
* ALL RIGHTS RESERVED.

THIS SOFTWARE IS FURNISHED UNDER A LICENSE AND MAY BE USED AND COPIED
ONLY IN ACCORDANCE WITH THE TERMS OF SUCH LICENSE AND WITH THE
INCLUSION OF THE ABOVE COPYRIGHT NOTICE. THIS SOFTWARE OR ANY OTHER
COPIES THEREOF MAY NOT BE PROVIDED OR OTHERWISE MADE AVAILABLE TO ANY
OTHER PERSON. NO TITLE TO AND OWNERSHIP OF THE SOFTWARE IS HEREBY
TRANSFERRED.

THE INFORMATION IN THIS SOFTWARE IS SUBJECT TO CHANGE WITHOUT NOTICE
AND SHOULD NOT BE CONSTRUED AS A COMMITMENT BY DIGITAL EQUIPMENT
CORPORATION.

DIGITAL ASSUMES NO RESPONSIBILITY FOR THE USE OR RELIABILITY OF ITS
SOFTWARE ON EQUIPMENT WHICH IS NOT SUPPLIED BY DIGITAL.

++
FACILITY:

LABIO demonstration system

ABSTRACT:

This module contains the I/O code for handling
an AD11-K. It is an example of a connect-to interrupt
routine. This module contains code to perform the following

The start I/O routine
The interrupt service routine
The cancel I/O routine

AUTHOR:

P. Programmer 15-Nov-79

--

.SBTTL DATA STRUCTURES
.PSECT LABIO_SECTION PIC,OVR,REL,GBL,SHR,NOEXE,RD,WRT,LONG

; The following data structures are also defined by a
; FORTRAN INCLUDE file. These definitions must agree.

; AD_BLOCK A/D Control Block

MAX_AD_CHANNEL = 16 ;Number of A/D channels
AD_BLOCK_SLOTS = 16 ;number of entries in one block
AD_BLOCK_SIZE = MAX_AD_CHANNEL*AD_BLOCK_SLOTS

;AD_BLOCK offsets (long words)

AD_STATUS	= 0	: STATUS (Unknown, inactive, or active)
ACTIVE	L= 2	: ACTIVE
INACTIVE	L = 1	: INACTIVE
PID	= 4	: PID of connected process
TICS_SAMPLE	= 8	: Rate in tics/sample
BUFFER_SIZE	= 12	: User specified buffer size
BUFFER_COUNT	= 16	: User specified buffer count
BUFFER_ACQ	= 20	: Number of buffers acquired
VALID_BUF_IND	= 24	: Index of current valid data buffer
VALID_BUF_COUNT	= 28	: Number of data points in last buffer
CUR_BUF_IND	= 32	: Index to current acq. buffer
CUR_BUF_COUNT	= 36	: Number of data points in last buffer
TICS_REMAINING	= 40	: Tics remaining to next sample
CUR_ACQ_OFF	= 44	: Offset to acq point
AD_BLOCK_END	= 64	: Offset to end of a block

AD_BLOCK: .BLKL AD_BLOCK_SIZE

; DATA_BUFFER Data buffers for LABIO

MAX_BUF_COUNT = 2 ;Number buffers/channel
MAX_BUF_SIZE = 512 ;Maximum buffer size (WORDS)

BUFFER_END = MAX_BUF_COUNT*MAX_BUF_SIZE*2 ; Size of one set of buffers

DATA_BUF_SIZE = MAX_AD_CHANNEL*MAX_BUF_SIZE*MAX_BUF_COUNT
DATA_BUFFER: .BLK# DATA_BUF_SIZE

DATA_BUFFER_OFF = DATA_BUFFER-AD_BLOCK ;Offset to data buffer from
;beginning of data structure

; CONNECT_BLOCK Process Connect control block

MAX_PID = 16 ;Max number of processes connected

CONNECT_SIZE = MAX_PID*2

CONNECT_BLOCK: .BLKL CONNECT_SIZE

LABIOCIN.MAR;1

16-SEP-1984 17:04:23.64 M⁹ Page 3

SC

.SBTTL I/O DEVICES

;This section defines the constants associated with the KW11-K clock
;and the AD11-K A/D converter

;KW11-K Clock
;CSR bit assignments
KW11\$M_GO = ^01 :GO bit
KW11\$M_RATE = ^02 :Rate = bits 2-4
KW11\$M_INTENB = ^0100 :Interrupt enable
KW11\$M_READY = ^0200 :Ready bit
KW11\$M_REPINT = ^0400 :repeated interuupts

KW11_CSR_CONS = KW11\$M_REPINT!KW11\$M_INTENB!<1*KW11\$M_RATE>
:Repeated interrupts,interrupt enable
:Rate = 1 MHz
KW11_PRESET = 1000. :Preset => Interrupt rate of 1 KHz

KW11_A_BUFFER = ^02 :Offset to clock A preset buffer
KW11_A_COUNTER = ^024 :Offset to clock A counter

;AD11-K A/D converter

AD11_OFFSET = -4 : Offset to the AD11 from thr KW11 clock CSR.
AD11_BUF = 2 : AD11 buffer offset fron AD11 CSR

AD11_GO = 1 : Go bit
AD11_MUX_INCR = ^0400 : Mux incr bit
AD11_CSR_CONS = AD11_GO : Initial CSR value

;Limit for stopping ISR loop
AD11_LOOP_LIMIT = AD11_MUX_INCR*<MAX_AD_CHANNEL-1>!AD11_CSR_CONS

\$IDBDEF : Definition for I/O drivers
\$UCBDEF : Data structurs
\$IODEF : I/O function codes
\$CINDEF : Connect-to-interrupt
\$CRBDEF : CRB stuff
\$VECDEF : more

.SBTTL LABIO_CIN_START, Start I/O routine

++
LABIO_CIN_START - Starts the KW11-K

Functional description:

This routine starts the KW11-K

Rate = 1 KHz

Repeated interrupt

Inputs:

0(R2) - arg count of 4

4(R2) - Address of the process buffer

8(R2) - Address of the IRP (I/O request packet)

12(R2) - Address of the device's CSR

16(R2) - Address of the UCB (Unit control block)

Outputs:

none

The routine must preserve all registers except R0-R2 and R4.

-- .PSECT LABIO_CIN

LABIO_CIN_START::

MOVL 12(R2),R0 ; Get address of the KW11 CSR
CLRW (R0) ; Clear the clock

MNEGW #KW11_PRESET,- ; Preset count buffer

MOVW #KW11_CSR_CONS+KW11\$M_GO,(R0) ; Set the bits for
; Repeated interrupt
; Interrupt Enable
; GO!

MOVW #SSS_NORMAL,R0 ; Load a success code into R0.
RSB ; Return

.SBTTL LABIO_CIN_INTERRUPT, Interrupt service routine

++
LABIO_CIN_INTERRUPT
Functional description:

Inputs:

0(R2) - arg count of 5
4(R2) - Address of the process buffer
8(R2) - Address of the AST parameter
12(R2) - Address of the device's CSR
16(R2) - Address of the IDB (interrupt dispatch block)
20(R2) - Address of the UCB (Unit control block)

Outputs:

Sets those bits in the AST parameter for those channels who had a buffer filled

The routine must preserve all registers except R0-R4

CIN_BUF_ADD = 4 ;Address of CIN buffer
AST_PARM = 8 ;Offset to AST parameter address
CIN_CSR_ADD = 12 ;Address of CSR

LABIO_CIN_INT::
PUSHR #^M<R5,R6> ;Service device interrupt, save R5,R6
MOVL CIN_CSR_ADD(R2),R4 ;Address of the KW11 CSR
MOVL CIN_BUF_ADD(R2),R5 ;Address of AD_BLOCK, control block
MOVAL DATA_BUFFER_OFF(R5),R1 ;for each A/D Channel
MOVAL AD11_OFFSET(R4),R4 ;Data Buffers
MOVL AD11_CSR_CONS,R6 ;Address of the AD11 CSR
CLRL #AST_PARM(R2) ;AD11 CSR bits, GO bit on
CLRL R3 ;Zero the AST parameter

AD_LOOP:
CMPL (R5),S^#ACTIVE_L ;Is this channel active?
BLSS AD_LOOP_NEXT ;No, try next channel

SOBGTR TICS_REMAINING(R5),AD_LOOP_NEXT
;Decr the timer for this channel
;Br if no conversion required

MOVW R6,(R4) ;Start conversion, while that's going on
.IF DF KW_HIST ;Time histogram, stored in data buffer
MOVZWL K011_A COUNTER-AD11_OFFSET(R4),R0 ;Get current clock contents
ADDW #KW1T_PRESET,R0 ;Calc time from intgerrupt
INCW (R1)[R0] ;Add one to that time bin
.ENDC

: While the A/D is converting, the tic counter for this channel,
: get the offset to the data pointer, and update it. Take appropriate
: action if we have buffer overflow.

```

MOVL TICS_SAMPLE(R5),- ;Reset timer for this channel
MOVL CUR_ACQ_OFF(R5),R0 ;Get index to next data point
INCL CUR_ACQ_OFF(R5) ;Advance it
AOBLSS BUFFER_SIZE(R5),- ;Update current data count
CUR_BUF_COUNT(R5),- ;Br if no buffer overflow
AD_LOOP_DATA

```

:Buffer overflowed, reset data pointer, reset buffer pointer
:increment acquired buffer count, termimate channel I/O if done

```

MOVL CUR_BUF_IND(R5),- ;Valid data buf available for user
MOVL CUR_BUF_COUNT(R5),- ;Number of points in buffer
MULL3 CUR_BUF_IND(R5),- ;Offset to next data point
#MAX_BUF_SIZE,-
CUR_ACQ_OFF(R5)
CLRL CUR_BUF_COUNT(R5) ;Reset data count
AOBLEQ #MAX_BUF_COUNT,- ;Next buffer index
CUR_BUF_IND(R5),1$ ;Wrap-around, reset buffer index
CLRL CUR_ACQ_OFF(R5) ;And buffer offset
1$: INSV #1,R3,#T,@AST_PARM(R2) ;Set bit in AST parameter word
AOBLSS BUFFER_COUNT(R5),- ;Incr buffer count
BUFFER_ACQ(R5),2$ ;Done with all buffers?
TSTL BUFFER_COUNT(R5) ;If original count was zero
BEQL 2$ ;Don't stop
MOVL #INACTIVE_L,(R5) ;Deactivate channel
2$:

```

: Now, get the data point and store it in the buffer.

```

AD_LOOP_DATA:
1$: TSTB (R4) ;Wait for A/D conversion
BGEQ 1$ ;
.IF NDF KW_HIST ;Time histogram don't store actual data
MOVW ADT1_BUF(R4),(R1)[R0] ;store data point in buffer.
.ENDC

```

:All done with this channel, setup for the next

```

AD_LOOP_NEXT:
ADDL #AD_BLOCK_END,R5 ;Next channel block
ADDL #BUFFER_END,R1 ;Next buffer
ADDW #AD11_MUX_INCR,R6 ;Incr A/D MUX
AOBLSS S^#MAX_AD_CHANNEL,R3,- ;Next channel
AD_LOOP ;Br if not done

```

```

:Exit routine - If any buffer overflowed, queue an AST
MOVL @AST_PARM(R2),R0 ;If any bit in the AST parameter
BEQL 1$ ;is set we must queue an AST
MOVL #1,R0 ;1 means queue the AST, 0 means don't
1$: POPR #^M<R5,R6> ;Restore R5,R6
RSB

```

.SBTTL LABIO_CIN_CANCEL, Cancel I/O routine

;+
; LABIO_CIN_CANCEL, Cancels an I/O operation in progress

Functional description:

This routine turns off the KW11-K

Inputs:
R5 - Addr of the UCB

Outputs:

The routine must preserve all registers except R0-R3.

;--

LABIO_CIN_CNCL::

MOVL	UCB\$L_CRB(R5),R0	: Get Address of the CRB
MOVL	CRB\$L_INTD+VE(\$L_IDB(R0),R0	:Address of the IDB
MOVL	IDB\$L_CSR(R0),R0	: Get addr of KW11
CLRW	(R0)	: Turn of the KW11
MOVW	#SSS_NORMAL,R0	: And return
RSB		

.SBTTL LABIO_CIN-END, End of Module
++
Label that marks the end of the module
--
LABIO_CIN-END: ; Last location in module

```
.SBTTL AD_CIN_SETUP Set-up routine for LABIO connect-to-interrupt
```

```
;+  
; This routine issues the QIO to connect to the AD11/KW11 interrupts.  
; It takes care of the internals associated with the connect-to-interrupt  
; QIO. Input parameters the VMS channel and the AST service routine address.  
; The connect-to-interrupt QIO condition code is returned.
```

```
.PSECT AD_CIN_SETUP
```

```
AD_CIN_SETUP::  
    .WORD 0  
    MOVL 8(AP),USER_AST ;Get the user AST routine addr  
AD_CIN_QIO:  
    $QIO_S CHAN=04(AP),- ;Channel  
    FUNC=#IOS_CININTWRITE,- ;Allow writing to the data buffer  
    IOSB=AD_CIN_IOSB,- ;I/O status Block  
    P1=AD_CIN_BUF_DESC,- ;Buffer descriptor  
    P2=#AD_CIN_ENTRY,- ;Entry list  
    P3=#AD_CIN_MASK,- ;Status bits,etc  
    P4=#AD_CIN_AST,- ;AST service routine  
    P6=#10 ;preallocate some AST control blocks  
    RET ;Return to caller
```

```
AD_CIN_BUF_DESC:  
    .LONG LABIO_CIN_END-AD_BLOCK ;Buffer descriptor for CIN  
    .LONG AD_BLOCK ;Size of buffer and CIN handler  
    .LONG LABIO_CIN_HANDLER ;Address of buffer
```

```
AD_CIN_ENTRY:  
    .LONG 0 ;No init code  
    .LONG LABIO_CIN_START-AD_BLOCK ;Start code  
    .LONG LABIO_CIN_INT-AD_BLOCK ;Interrupt service routine  
    .LONG LABIO_CIN_CANCEL-AD_BLOCK ;I/O cancel routine
```

```
AD_CIN_IOSB:  
    .LONG 0,0 ; I/O Status Block
```

```
; Control mask
```

```
AD_CIN_MASK = CIN$M_REPEAT!CIN$M_START!CIN$M_ISR!CIN$M_CANCEL
```

```
; AD_CIN_AST
```

```
; This AST routine calls the user AST routine. The user routine  
; can not be called directly because the AST parameter itself  
; not its address is returned via the connect-to-interrupt routine.  
; This routine simply calls the user routine with the ADDRESS of  
; the AST parameter.
```

```
AD_CIN_AST::  
    .WORD 0  
    PUSHAL 4(AP) ;Get the AST parameter addr  
    CALLS #1,USER_AST ;Call the USER routine  
    RET
```

USER_AST:

```
.LONG 0  
.END
```

;Addr of the user AST routine

0157 AH-BT13A-SE
VAX/VMS V4.0

DIGITAL EQUIPMENT CORPORATION
CONFIDENTIAL AND PROPRIETARY

LPMULT
B32

DRMAST
MAR

XADRIVER
MAR

TORIVER
MAR

USSTEST
MAR

GBLSECUFO
MAR

USSDISP
MAR

DOO ERAPAT
MAR

LBRMAC
MAR

XADRIVER
MAR

WORKO
LIS

LABIUCIN
MAR

SCRF
MAR

DRSLU
MAR

DTE DF00
MAR

EXAMPLES